

KOALA NUTRITIONAL WARMER™



Designed by a Former NICU Nurse

The Koala is an easy-to-use UL listed nutritional warmer with an integrated temperature feedback loop. The feedback loop is operated by custom algorithms that prevents damage to the properties of human breast milk and fortifiers. The Koala's intuitive user interface guides the operator through the three modes of warming and thawing "Fridge to Feed", "Frozen to Feed" and "Frozen to Thaw".



| Catalog Number | Description | Unit of Measure |
|----------------|--------------------------|-----------------|
| APG-KNW-1 | Koala Nutritional Warmer | 1/ea |
| APG-K-POUCH | Koala Pouch | 300/ca |



The Kentec Koala™ Nutritional Warmer

ADVANTAGES

The CDC recommends thawing milk the following ways:

- In the refrigerator overnight
- In container of warm or lukewarm water
- Under lukewarm running water

“Water will transfer heat 24.17 times faster than air.”¹

Koala Pouches are delivered 300 bags per case. 6 boxes of 50 bags each.

The red LED numbers tells you how long that amount of nutrition took to reach proper temperature.

The green LED time tells you how much time is left to use the nutrition (counts down from 30 minutes).

When the green LED number (Ready: 30) appears along with the green LED light bar the nutrition is ready, and the 30-minute countdown has started.

The heater feedback loop and custom algorithms prevent nutrition from being overheated.

The constant vibration of the container throughout the entire warming cycle provides evenly heated nutrition.

UL Listed, 61010-1



¹forum.heatinghelp.com/discussion/154525/water-vs-air-as-a-heat-transfer-medium Harvey Ramer

Warmer Comparison

| | Kentec Koala | Ameda* | Medela** |
|---|--------------|--------|----------|
| Uses vibration to mix nutrition for even heating | ● | ● | |
| Warmer incorporates custom algorithms | ● | ● | ● |
| Warm from “Fridge to Feed” temperature | ● | ● | ● |
| Warm from “Frozen to Thaw” temperature | ● | ● | ● |
| Warm from “Frozen to Feed” temperature | ● | ● | ● |
| Easy to clean | ● | ● | ● |
| Accommodates baby bottles up to 250ml and syringes from 1 to 100ml | ● | ● | |
| Compatible with most brands of breast milk storage bags | ● | ● | |
| Warmer uses a heater with sensor feedback loop | ● | ● | |
| Dual pouch bag design | ● | ● | |
| The extended fold over flap prevents water from splashing into the dedicated syringe or bottle area while filling | ● | | |
| The warmer interface guides user through an easy 3-step process | ● | | |
| Red LED lights indicate the nutrition is heating | ● | | |
| Green LED numbers and light bar indicate the nutrition is ready and a 30- minute count down has started | ● | | |

*Ameda is a registered trademark of Ameda, Inc. 485 Half Day Rd Buffalogue, FL 60084 USA

**Medela is a registered trademark of Medela, Inc. 1101 Corporate Dr. McHenry, IL 60050 USA