



Medfusion® 3500 Syringe Pump and PharmGuard® Medication Safety Software

Features Ameritus® enteral syringe recognition

Reduce the risk of IV/enteral misconnections

By partnering with Kentec Medical, Smiths Medical ensures you'll have the best of both worlds - a smarter infusion system and safer enteral delivery.



Ameritus® Oral Enteral



Ameritus® ENFit® Enteral

Medfusion® 3500 syringe pumps with PharmGuard® medication safety software now recognize the following Ameritus® syringe sizes:

- 5mL
- 10mL
- 20mL
- 30mL
- 60mL

- Orange color faceplate for easy enteral identification*
- Unique tips will not connect to standard IV Luers

Adding enteral syringe data to PharmGuard® medication safety software** allows you to:

- Comply with Joint Commission recommendations for enteral feedings¹
- Differentiate your patients' IV and enteral lines²
- Meet your facility's enteral delivery needs

Ameritus® is a federally registered trademark of Kentec Medical, Inc. ENFit® is a federally registered trademark of GEDSA, used with permission.

* Available through third party supplier

**Requires software maintenance for all Medfusion® 3500 and Medfusion® 4000 pumps in the facility.

1. Guenter, Peggy, PhD., R.N., C.N.S.N., et al. Enteral Feeding Misconnections: A Consortium Position Statement. The Joint Commission Journal on Quality and Patient Safety. May 2008; Volume 34 Number 5

2. Abstract: Enteral feeding set connectors and adapters. Association for the Advancement of Medical Instrumentation. ANSI/AAMI ID54: 1996/© 2005

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN CANADA AND OTHER COUNTRIES

PSS - SMF3500 Ameritus rev New

Contact Us: (800) 825-5996 | Email: sales@kentecmedical.com

Kentec Medical Inc. | Our Products Provide a Lifetime of Care

(949) 863-0810 | www.kentecmedical.com

